



W. W. Hansen Experimental Physics Laboratory
STANFORD UNIVERSITY
STANFORD, CALIFORNIA 94305 - 4085

Gravity Probe B Relativity Mission

PROCEDURE FOR
Cryo Focal Position Measurement of Science Telescope
For Cryo Focal Test

GP-B P0252 Rev -

June 18, 1997

Prepared by: Suwen Wang
Engineer

Date

Approved by: John Lipa
Manager, Telescope Development

Date

Approved by: B. Taller
Quality Assurance

Date

Approved by: J. Turneure
Hardware Manager

Date

GP-B Procedure Document 252

Cryo Focal Position Measurement for Science Telescope For Cryo Focal Test

RE: Suwen Wang
DATE: June 18, 1997
ESTIMATED DURATION: 1 week

Objective:

To determine the focal position of Science Telescope at 4.2 K.

Requirements:

- Procedure to be performed by certified personnel only.
Certified personnel for the test:
Suwen Wang
- Procedures as described in Pdoc P0251 complete.

Calibrations:

- The artificial star is calibrated for its collimation using the self contained knife edge device.
- Other measurements are relative and require no calibration.

Precautions:

- Science Telescope is made of fused quartz, a fragile material. Extra precaution is to be taken during all handling to prevent damaging the Telescope.

Conditions for the test:

- Telescope under test:
Dwg No: 25091-204 Rev - _____.
Telescope Serial No. _____.
- Test to be performed in artificial star #2 lab.

1. Procedures for telescope probe/test dewar integration:

Lift the probe up by 1/2" with the hoist and move away the probe cart.

Assemble the copper tube to the outside of the vacuum can.

Attach heater and thermometer wires.

Use masking tape to protect the heater and thermometer wire pin connectors.

Lower the probe down into the dewar carefully to ensure the centration of the probe o.d.

with respect to the dewar i.d.

Tie down ten bolts which holds the probe onto the dewar with 6 in-lb torque.

Remove the lifting plate from top of the probe by unscrew ten screws.

Signed **Date**

- Procedures for attaching the probe/dewar assembly to the artificial star: _____
- Connect the cables for the readout electronics. _____
- Use a flash light to direct light into the telescope and look for signal changes in the detector readout to verify the functionality of the detectors. _____
- Pump the suspension tires to about 50 psi. _____
- Manipulate the dewar/probe assembly to beneath the star. _____
- Align the screw holes on top of the probe to those on the star. _____
- Deflate the suspension till the gap of the star and telescope is about 1/4". _____
- Place 10 1/4-20 screws from below top of the probe. _____
- Further deflate the suspension till the gap is almost zero. _____
- Tighten all the screws. _____
- Pump the suspension back to 55 psi. Now the dewar and telescope probe assembly should be suspended from the floor. _____
- Align the star beam such that the focused beam would fall on to two of the 4 detectors of the quad across the direction of the intended scan. _____
- _____
- _____
- _____
- _____

4. Measurement of the focal position

4.1 First, the focal position is measured at room temperature by changing the collimation of the star such that it is equivalent to the detector being +/- 0.15" from the telescope focal point along the optical axis and monitor the slope in the mid portion of the scan curve. Such measurement is achieved by adjusting the fiber tip location of the artificial star light source and running my software FocusScan.

The data were recorded in files named scan_Dir#_date
 where:

Dir can be either x or y for the scan direction

is the serial number of the scan of the day

date is the date in the format of m/d/y

The data format of the file is in column vectors. They are ordered in the following manner:

x(arc sec) y(arc sec) chan1(V) chan2(V) chan3(V) chan4(V)

y column is always 0 in this test.

Enter the file names in Table 1.

4.2 The telescope is then cooled to 4.2K.

5.6 Procedure 5 complete.

Signed: _____

Date: _____